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Warenzeichen / Art.-Nr.

Alfa Nitril-Einweghandschuh

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766500M Size 8
766500L Size 9
76650XL Size 10

Description

Thin nitrile disposable glove with beaded cuff. The glove is powder free and ambidextrous.

General information

Before use, it should be tested/ensured that the glove has the appropriate size to achieve the best possible comfort and safety at work.

The service life cannot be determined and depends on the scope of application and the extent to which the user makes sure that the glove is suitable for the intended use. No allergy-causing substances are known.

Category

The glove is certified in category III in compliance with Medical Device Directive 93/42/EEC and Personal Protective Equipment Regulation (EU) 2016/425.

The glove has been tested in accordance with the standards EN 420:2003+A1:2009 (general requirements), EN 455 (medical use) and EN ISO 374-1:2016, EN 374-4:2013, EN ISO 374-5:2016 (chemicals and microorganisms).

Notified body responsible for EU Type Examination and Module C2 ongoing conformity assessment: SATRA Technology Europe Ltd, Bracetown Business Park, Clonee, Dublin, D15 YN2P, Ireland (Notified Body No. 2777).

EN 374 Test Data / Type B

EN ISO 374-1:2016

	Chemical	Time*	Level	Degradation %
A	Methanol	N/A	N/A	N/A
B	Acetone	N/A	N/A	N/A
C	Acetonitrile	N/A	N/A	N/A
D	Dichloromethane	N/A	N/A	N/A
E	Carbon disulphide	N/A	N/A	N/A
F	Toluene	N/A	N/A	N/A
G	Diethylamine	N/A	N/A	N/A
H	Tetrahydrofurane	N/A	N/A	N/A
I	Ethyl acetate	N/A	N/A	N/A
J	N-heptane	N/A	N/A	N/A
K	Sodium hydroxide 40%	>480	6	-42.9
L	Sulphuric acid 96%	N/A	N/A	N/A
M	Nitric acid 65 %	<10	0	97.6
N	Acetic acid 99 %	<10	0	93.9
O	Ammonium hydroxide 25%	>10	1	-52.0
P	Hydrogen peroxide 30%	>30	2	22.8
S	Hydrofluoric acid 40 %	N/A	N/A	N/A
T	Formaldehyde 37 %	>60	3	5.0
	**4% Chlorhexidine Digluconate	>480	6	19.0
	10-13% Sodium Hypochlorite	>480	6	14.7
	50% Sulphuric acid	>480	6	-20.5
	10% Acetic acid	>120	4	66.7
	5% Ethidium Bromide	>480	6	3.4
	50% Glutaraldehyde	>480	6	27.4

0.1% Phenol	>480	6	33.8
1.5% Methanol in water	>480	6	21.9
70% Isopropanol	<10	0	62.2
35% Ethanol	<10	0	38.8
3% Povidone-iodine	>480	6	33.7
10% Sodium Percarbonate	>480	6	15.4

*)Breakthrough time (minutes)

***)Permeation rate: 7 µg / cm² / min

Performance Level	1	2	3	4	5	6
Minimum break through time (mins)*	>10	> 30	> 60	> 120	> 240	> 480

*) Glove performance quoted is based on laboratory data and may not reflect the actual duration of protection in the workplace due to other factors influencing the performance such as temperature, abrasion, degradation etc.

EN374-4:2013

Degradation levels indicate the change in puncture resistance of the gloves after exposure to the challenge chemical.

EN ISO 374-5:2016

Resistance to Bacteria and Fungi = Pass

Resistance to Virus = Pass

Migration test

(EN 1186-1:2002)

Application / qualities

Thin and flexible seamless glove for great dexterity. The glove is waterproof and suitable for handling of many different chemicals with short contact time. If the glove has been in contact with substances it must be discarded when the breakthrough time is reached. Primarily used in food production, mounting, product control, laboratories, chemical industry and for general cleaning.

Treatment / storage

The glove is for single use and should not be cleaned or re-used. Best stored in dark and cool surroundings in the original packaging.

Warning

- This information does not reflect the actual duration of protection in the workplace and the differentiation between mixtures and pure chemicals.
- The chemical resistance has been assessed under laboratory conditions from samples taken from the palm only and relates only to the chemical tested. It can be different if the chemical is used in a mixture.
- It is recommended to check that the gloves are suitable for the intended use because the conditions at the workplace may differ from the type test depending on temperature, abrasion and degradation.
- When used, protective glove may provide less resistance to the dangerous chemical due to changes in physical properties. Movements, snagging, rubbing, degradation caused by chemical contact etc. may reduce the actual use time significantly. For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistant gloves.
- Before usage, inspect the gloves for any defect or imperfections.
- The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimens.

Packaging

100 gloves per dispenser.

10 dispensers in a carton made of recyclable cardboard.

Alfa GmbH

Ferdinand-Porsche-Straße 10

73479 Ellwangen / Germany

Tel.: +49 (0) 7961-57 99 0

alfa-direct.com

CE 2777

KPT

VIRUS