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## Warenzeichen / Art.-Nr.

### Alfa Diamond Touch

9430050M Größe M/8

9430050L Größe L/9

9430045XL Größe XL/10

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9430045XL Size XL/10

### Description

Thin nitrile disposable glove with beaded cuff and diamond structure. The glove is powder free and ambidextrous.

### General information

Before use, it should be tested/ensured that the glove has the appropriate size to achieve the best possible comfort and safety at work.

The service life cannot be determined and depends on the scope of application and the extent to which the user makes sure that the glove is suitable for the intended use.

### Category

The glove is certified in category III in compliance with PSA-regulation (EU) 2016/425.

The glove has been tested in accordance with the standards EN 420:2003 + A1:2009 (general requirements) and EN 374-1:2016, EN 374-2:2003 (AQL 1,5 (Level 2)), EN 374-5:2016 (chemicals and micro organisms). The glove is also EN 455-1-2-3 tested.

### EN ISO 374-1:2016/Type B

|          | Chemical              | Level |
|----------|-----------------------|-------|
| <b>K</b> | Sodium Hydroxide 40 % | 6     |
| <b>A</b> | Methanol              | 0     |
| <b>G</b> | Diethylamin           | 0     |
| <b>O</b> | Chlorhexidine 4 %     | 6     |
| <b>T</b> | Formaldehyde 37 %     | 6     |
| <b>T</b> | Glutaraldehyde 1%     | 6     |
| <b>T</b> | Glutaraldehyde 4 %    | 6     |
| <b>P</b> | Hydrogen peroxide 3 % | 6     |
| <b>T</b> | Formalin 10 %         | 6     |
| <b>J</b> | n-Heptane             | 2     |
| <b>I</b> | Ethyl Acetate         | 0     |
| <b>F</b> | Toluene               | 0     |
| <b>D</b> | Dichloromethane       | 0     |
| <b>E</b> | Carbon Disulphide     | 0     |

| Performance Level                 | 1  | 2  | 3  | 4   | 5   | 6   |
|-----------------------------------|----|----|----|-----|-----|-----|
| Minimum break through time (mins) | 10 | 30 | 60 | 120 | 240 | 480 |

Glove performance quoted is based on laboratory data and may not reflect the actual duration of protection in the workplace due to other factors influencing the performance such as temperature, abrasion, degradation etc.

### EN ISO 374-5:2016

Resistance to Bacteria and Fungi = Pass

### Migration test

(EN 1186-1:2002)

### Application / qualities

Thin and flexible seamless glove for great dexterity.

The glove is waterproof and suitable for handling of many different chemicals with short contact time.

If the glove has been in contact with substances it must be discarded when the breakthrough time is reached.

Primarily used in mounting.

### Treatment/storage

The glove is for single use and should not be cleaned or re-used. Best stored in dark and cool surroundings in the original packaging.

### Warning

- This information does not reflect the actual duration of protection in the workplace and the differentiation between mixtures and pure chemicals.
- The chemical resistance has been assessed under laboratory conditions from samples taken from the palm only and relates only to the chemical tested. It can be different if the chemical is used in a mixture.
- It is recommended to check that the gloves are suitable for the intended use because the conditions at the workplace may differ from the type test depending on temperature, abrasion and degradation.
- When used, protective glove may provide less resistance to the dangerous chemical due to changes in physical properties. Movements, snagging, rubbing, degradation caused by chemical contact etc. may reduce the actual use time significantly. For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistant gloves.
- Before usage, inspect the gloves for any defect or imperfections.

### Packaging

100/90 gloves per dispenser.

10 dispensers in a carton


### Alfa GmbH

Ferdinand-Porsche-Straße 10

73479 Ellwangen / Germany

Tel.: +49 (0) 7961-57 99 0

alfa-direct.com

CE 0321

EN ISO 374-1:2016/Type B

EN ISO 374-5:2016

KOTP